

Certificate of Compliance



Part **Piston Pump CP Xtend™**
Part Number **12860 (6-channels)
12870 (7-channels)
12880 (8-channels)**
Year of Construction **from 2013-08**

This document certifies that the above mentioned apparatus was assembled and successfully tested at the SOTAX AG facility in Aesch, Switzerland.

- DQ = Design qualification
- Visual check of system
- Electrical safety check
- Mechanical testing
- Calibration
- System testing
- Compatibility with the valid USP and Ph. Eur. requirements

We confirm that all parts in contact with the dissolution media are made of inert material such as:

Tubes for medium	PTFE (optional PEEK)
Pump heads, connecting blocks	PVDF
Piston, cylinder	Ceramics
Gaskets	VITON
Pump head gaskets	UHMWPE
Connection fittings	FPM

Test data for this apparatus is available for inspection on request.

Aesch, 2019-07-26

Bruno Balli
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland