

Certificate of Compliance



Part DT50
Part Number 14770-01 (Master Station)
14770-02 (Client Station)
14780-01 (MediaPrep Station)

This document certifies that the above mentioned apparatus was assembled and successfully tested at the SOTAX AG facility in Aesch, Switzerland.

- DQ = Design qualification
- Visual check of system
- Electrical safety check
- Mechanical testing
- Calibration
- System testing
- Compatibility with the valid USP and Ph. Eur. requirements

We confirm that all parts in contact with liquid are made of inert material such as:

Basket holder and plate	Polymethyl methacrylate (PMMA), Polyoxymethylene (POM), Polyvinyl chloride (PVC), Polyamide (PA), Stainless steel
Mesh	Stainless steel (SB) / PE (SK)
Beaker	Borosilicate Glass
Metal plate for induction heater	Stainless steel

Test data for this apparatus is available for inspection on request.

Aesch, 2023-05-03

Jürg Suter
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland