



Part **Basket without Lock**  
 Part Number **02830**  
 Serial Number **N11017764 - N11018284**



# Test Certificate

|                   |   | <b>USP / Ph.Eur.<br/>Tolerances</b> |         |                |
|-------------------|---|-------------------------------------|---------|----------------|
| <b>Dimensions</b> | Basket                                  | Basket rim O.D. Ø                   | 25.0 mm | ± 3.0 mm       |
|                   |   | Height                              | 37.0 mm | ± 3.0 mm       |
|                   |   | Wire Ø                              | 0.25 mm | 0.22 – 0.31mm  |
|                   |   | Wire openings                       | 0.4 mm  | 0.36 - 0.44 mm |
|                   |   | Screen opening height               | 27.0 mm | ± 1.0 mm       |
|                   |   | Screen O.D. Ø                       | 22.2 mm | ± 1.0 mm       |
|                   |   | Clear opening Ø                     | 20.2 mm | ± 1.0 mm       |
| <b>Material</b>   | Stainless Steel (electrolytic polished) |                                     |         |                |

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Aesch, 2023-04-10

Michael Schmidlin  
 SOTAX Quality Management  
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland