



Part Suction Head with Stirring Shaft
 Part Number 07570
 Serial Number N15000544 - N15000868



Test Certificate

			USP / Ph.Eur. Tolerances
Dimensions	Shaft	Diameter	9.4 – 10.1 mm
Material	Suction head PVC		
	Shaft Stainless Steel (electrolytic polished)		
	Suction tube PTFE		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Aesch, 2019-03-06

Bruno Balli
 SOTAX Quality Management
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland