



Part TTS-Rotor with Stirring Shaft  
(Length 50.8 mm)

Part Number 07663

Serial Number N30000001 - N30000030



# Test Certificate

			USP / Ph.Eur. Tolerances
<b>Dimensions</b>	Length	50.8 mm	± 0.1 mm
	Outer Ø	44.5 mm	- 0.1 mm
	Wall thickness	0.75 mm	± 0.1 mm
	4 holes Ø	11.1 mm	± 0.1 mm
	Bolt circle (equally spaced)	24.5 mm	± 0.1 mm
	Angle to surface	63.4 °	± 0.5 °
	Shaft Ø	9.4 – 10.1 mm	- -
<b>Material</b>	Stainless Steel (electrolytic polished)		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Aesch, 2019-01-07

Jürg Suter  
SOTAX Quality Management  
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland