


# Test Certificate

Part: **Disintegration Glass Tube**

Number: 08825-1

For SOTAX DT Tablet Disintegrators



Dimensions:	a	Glass strength	2.5 mm	Ph.Eur. Tolerances
	d	Inside Ø	33.0 mm	± 0.5 mm
	l	Length	77.5 mm	± 0.5 mm
				± 2.5 mm

Material: Glass

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with Ph. Eur. 2.9.1.-2 Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management  
Bruno Balli  
2018-07-23

  
Signature