




Part 1 Liter Vessel
 Part Number 14143-11
 Serial Number P35 001959 – P35 002101



Test Certificate

		USP / Ph.Eur. Tolerances		
Dimensions	Vessel	Inside Ø a	102.9 to 105.1 mm	98 – 106 mm
		Height b	200 mm	160 – 210 mm
		With hemispherical bottom		
Material	Borosilicate glass			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Navi Mumbai, 2023-02-16

16-02-2023 14:30:18

Moreshwar Pulekar
 SOTAX Quality Management
 C-64 | First Floor | Suchi Business Park | TTC Industrial area | MIDC Turbhe | Navi Mumbai 400705 | India