




Part Mini Vessel  
Part Number 14535-01  
Serial Number 97842 - 97848



		Tolerances	
Dimensions	Vessel	Inside Ø a	40 mm ± 1 mm
		Height b	203.8 mm ± 0.5 mm
		With hemispherical bottom	
Material	Borosilicate glass		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed.

Aesch, 2024-03-07

Michael Schmidlin  
SOTAX Quality Management  
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland

# Test Certificate