




Part Automated Disintegration Disk Set (6 x Disk)
 Part Number 23441-01



Test Certificate

USP / Ph.Eur. Tolerances

Specific gravity	ρ 1.18 – 1.20 kg/dm ³	--
Material	PMMA	

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Aesch, 2021-02-25

Bruno Balli
 SOTAX Quality Management
 SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland