


Test Certificate

Part: **Beaker**
Number: 03061-2
Manufacturer: SOTAX
Serial Number: 54755 - 54764

For SOTAX DT Tablet Disintegrators



Dimensions:

Inside Ø
Height

USP/Ph.Eur.
Tolerances
97 - 115 mm
138 - 160mm

With straight bottom

Material:

Borosilicate glass

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management
Bruno Balli
2013-12-09


Signature