


# Test Certificate

Part: **Beaker**  
Number: 03061-2  
Serial Number: 70887 - 70936

For SOTAX DT Tablet Disintegrators



Dimensions:

Inside Ø  
Height

USP/Ph.Eur.  
Tolerances  
97 - 115 mm  
138 - 160mm

With straight bottom

Material:

Borosilicate glass

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management  
Bruno Balli  
2017-06-07

  
Signature