

Certificate of Compliance



Part AT Xtend™
Part Number 15160, 15170, 15180 (1 Liter Units)
15260, 15270, 15280 (2 Liter Units)
Year of Construction from 2014-09

This document certifies that the above mentioned apparatus was assembled and successfully tested at the SOTAX AG facility in Aesch, Switzerland.

- DQ = Design qualification
- Visual check of system
- Electrical safety check
- Mechanical testing
- Calibration
- System testing
- Compatibility with the valid USP and Ph. Eur. requirements

We confirm that all parts in contact with the dissolution media are made of inert material such as:

Tubes for medium with active substance	PTFE (Teflon) PEEK (optional)
Test vessel	Borosilicate glass Polycarbonate (optional)
Suction head	PVDF
Seals	VITON (FPM)
Various	Stainless steel: A316L resp. 1.4401, 1.4404, 1.4435, 1.4460, 1.4462

Test data for this apparatus is available for inspection on request.

Aesch, 2019-07-30

Bruno Balli
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland