



Part Basket without Cover, 40 Mesh
 Part Number 02830
 Serial Number N11020355 - N11020426



Test Certificate

		USP / Ph.Eur. Tolerances		
Dimensions	Basket	Basket Lower Ring O.D. Ø	25.0 mm	± 3.0 mm
		Height	37.0 mm	± 3.0 mm
		Wire Ø	0.25 mm	0.22 – 0.31mm
		Wire openings	0.4 mm	0.36 - 0.44 mm
		Screen opening height	27.0 mm	± 1.0 mm
		Screen O.D. Ø	22.2 mm	± 1.0 mm
		Clear opening Ø	20.2 mm	± 1.0 mm
Material	Stainless Steel (electrolytic polished)			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements of a random sampling of the serial number range are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Westborough, 2024-04-11



11.04.2024 15:35:40

Jeff Schiavoni
 SOTAX Quality Management
 SOTAX Corp | 2400 Computer Dr | Westborough, MA 01581 | USA