



# Test Certificate

Part **Basket without Lock**  
Part Number **02830**  
Serial Number **N11003465 – N11003682**



		USP / Ph.Eur. Tolerances		
<b>Dimensions</b>	Basket	Basket rim O.D. Ø	25.0 mm	± 3.0 mm
		Height	37.0 mm	± 3.0 mm
		Wire Ø	0.25 mm	0.22 – 0.31mm
		Wire openings	0.4 mm	0.36 - 0.44 mm
		Screen opening height	27.0 mm	± 1.0 mm
		Screen O.D. Ø	22.2 mm	± 1.0 mm
		Clear opening Ø	20.2 mm	± 1.0 mm
<b>Material</b>	Stainless Steel (electrolytic polished)			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopoeia requirements (at date of approval, see below).

Aesch, 2018-12-12

Bruno Balli  
SOTAX Quality Management  
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland