



Test Certificate

Part Suction Head with Filter Insert (500 ml /900 ml)
Part Number 07760
Serial Number N19000000 – N19000089



			USP / Ph.Eur. Tolerances
Dimensions	Shaft	Diameter	9.4 – 10.1 mm
Material	Suction head PVC		
	Shaft Stainless Steel (electrolytic polished)		
	Suction tube PTFE		

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with USP and Ph. Eur. Pharmacopeia requirements (at date of approval, see below).

Aesch, 2020-01-31

Bruno Balli
SOTAX Quality Management
SOTAX AG | Nordring 1 | 4147 Aesch | Switzerland