


Test Certificate

Part: **Disintegration PMMA Tube**

Number: 08825-2

For SOTAX DT Tablet Disintegrators



				Ph.Eur. Tolerances
Dimensions:	a	Glass strength	2.5 mm	± 0.5 mm
	d	Inside Ø	33.0 mm	± 0.5 mm
	l	Length	77.5 mm	± 2.5 mm
Material:	PMMA			

The accuracy of the above mentioned part dimensions have been examined. All relevant measurements are within tolerance and the factory acceptance test has passed. The dimensions are in conformance with Ph. Eur. 2.9.1.-2 Pharmacopeia requirements (at time of delivery).

SOTAX Quality Management
Bruno Balli
2018-07-23


Signature